


Total IV Device Disinfection



SwabArt[®] Disinfecting Cap

SwabArt® Disinfecting Cap for swabable needle-free valves is a single-use disinfecting product that protects patients from infection.

The product is impregnated with a high liquid absorbent sponge saturated with 70% Isopropyl Alcohol (IPA).



- Provides active disinfection for up to 7 days if not removed
- SwabArt® is compatible with the EN ISO 80369-7-standard
- Compatible with all swabable needle-free valves
- All materials are sterile, and latex-free
- Strip design allows hanging on drip poles
- Supplied in boxes of 25 strips of 10 disinfecting caps
- Supplied in boxes of 250 individual disinfecting caps

Code	Product Description	Colour	Form	# of Caps	# of Strips/Box	Total Caps/Box
C100	SwabArt® Disinfecting Cap	Turquoise	Strip	10	25	250
C200	SwabArt® Disinfecting Cap	Turquoise	Individual Cap	1	-	250

“One CLABSI prevention method is disinfection of intravenous access points..... Use of a disinfectant cap on IV needleless connectors.....was associated with decreased CLABSI and costs.”

“Disinfectant cap use was associated with an estimated savings of almost \$300,000 per year in the hospital studied.”

Merrill, Katreena Collette et al. American Journal of Infection Control, 42, 12, 2014, p1274 - 1277

“...easy to use...”, “..With the alcohol-impregnated caps, it’s one less thing to think about.”, “Better Safety, Better Resource Utilization”

Wawrzyniak, Marcelina M., American Journal of Infection Control, 42, 6, 2014, P S15-S16