


SwabArt® Disinfecting Cap



Product Code
C100-C200

SwabArt® Disinfecting Cap User's Guide

The SwabArt® Disinfectant Cap is used to prevent the risk of infectious disease by providing continuous disinfection of the needle-free valve surface



SwabArt® Disinfecting Cap
Product Video



The cap may remain on the valve for seven (7) days if not removed



Do not use the cap if the seal/foil is damaged or detached



To be used only on needle-free valves



SwabArt® Disinfecting Cap is a single-use product do not re-use or re-sterilise