

Total IV Device Clarity




FlowArt[®] Needle-free Valves

FlowArt® Sterile Infusion Therapy Devices

KMEN.006_Rev00/102020

FlowArt® is a needle-free valve that has a **fully transparent clear housing** with an integrated flat silicone seal and internal fluid pathway, that protects the patient and nursing staff from exposures.



Product Code	Product Image	Product Description	Units per box	Priming Volume (ml)	Length (cm)	Flow Rate (ml/min)
A1010		FlowArt® Needle-free Valve	250	0,085	2,5	312

Manufacturing Facility

İkitelli OSB Mahallesi, 17. Cadde, No:17
Basaksehir, Istanbul - TURKEY

Phone : +90 212 494 27 27
Fax : +90 212 494 34 37

www.assetmedikal.com
info@assetmedikal.com

asset
Design the future