

NEUTRAL DISPLACEMENT

Needle-Free Connector



NeutrArt[®]
Sterile Infusion Therapy Devices

A1010N

NeutrArt® Needle-free Valve Neutral



Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min

AU1010N

NeutrArt® Needle-free Valve Neutral with Extension Line



Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

AUWB1010N


NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line



Priming Volume	Length	Flow Rate
0.46 ml	10 cm	90 ml/min

A2010N


NeutrArt® Double Lumen Needle-free Valve Neutral



Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010N

NeutrArt® Triple Lumen Needle-free Valve Neutral



Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

AU1010NWOR


NeutrArt® Needle-free Valve Neutral with Extension Line, Luer Slip



Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

A2010NWOR


NeutrArt® Double Lumen Needle-free Valve Neutral, Luer Slip



Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010NWOR

NeutrArt® Triple Lumen Needle-free Valve Neutral, Luer Slip



Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

ATU3010N

NeutrArt® Needle-free Valve with Three Way Stopcock, Neutral



Priming Volume	Length	Flow Rate
0.48 ml	16.5 cm	60 ml/min

A1010N-PH

NeutrArt® Needle-free Valve Neutral, Photoprotective



Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min

AU1010N-PH

NeutrArt® Needle-free Valve Neutral with Extension Line, Photoprotective



Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

AUWB1010N-PH


NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line, Photoprotective



Priming Volume	Length	Flow Rate
0.46 ml	10 cm	91 ml/min

A2010N-PH


NeutrArt® Double Lumen Needle-free Valve Neutral, Photoprotective



Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010N-PH

NeutrArt® Triple Lumen Needle-free Valve Neutral, Photoprotective



Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

AU1010NWOR-PH


NeutrArt® Needle-free Valve Neutral with Extension Line, Luer Slip, Photoprotective



Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

A2010NWOR-PH


NeutrArt® Double Lumen Needle-free Valve Neutral, Luer Slip, Photoprotective



Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010NWOR-PH

NeutrArt® Triple Lumen Needle-free Valve Neutral, Luer Slip, Photoprotective



Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min



Custom length sizes available upon request.

ASSET® MEDİKAL TASARIM SAN. VE TİC. A.Ş.
İkitelli OSB Mahallesi, 17. Cadde, No:17
Basaksehir / İSTANBUL TURKEY
T: +90 212 494 27 27 F: +90 212 494 34 37
info@assetmedikal.com - www.assetmedikal.com

asset
design the future

KMEN.009_Rev00/012021