

FLOWART® Y SITE NEEDLE-FREE VALVE

DESCRIPTION


FlowArt® Y Site Needle-free Valve is a connector with an integrated silicone seal valve.

FEATURES

- ✓ Microbial barrier swabable silicon stem with a split septum, flat seal.
- ✓ Clear transparent housing to allow clear view into the fluid pathway for enhanced flushing
- ✓ Compatible with TPN solutions, blood and oncology drugs.
- ✓ Suitable for high pressure infusions (up to 325 psi).
- ✓ Allows multiple accesses.
- ✓ ISO 80369-7 compatible.
- ✓ MRI and CT suitable.
- ✓ Shelf Life: 5 Years
- ✓ in compliance with MDD 93/42CEE. EC Certificate No: 1994-MD-18-524
- ✓ Available in either bulk / OEM and packed / sterile.
- ✓ Photoprotective (-PH) versions are in amber color.
- ✓ Latex-free, PVC-free, DEHP-free

Part	Material
Valve Body	Copolyester
Swabable septum	Silicone

DIAGRAM



PRODUCT DATA

Characteristic	Specification	
Product Code	BM04-39, BM04-PH-39	BM04 BM04-PH
Product Name	Needle-free Valve Port	Needle-free Valve Port
Priming Volume (mL)	0.49	0.50
Flow Rate (mL/min @ 1m water head pressure)	290	290
Flushing Volume (mL)	0.70	0.75
Length (cm)	3.5	3.5

*blue lower part

LABELLING AND PACKAGING

Each unit is labelled on its blister and packaging container.

Each container is labelled with the following traceability information:

- Quantity,
- Product description,
- Product date,
- Lot number

TDS.112_Rev03_26.04.2021

ASSET® MEDICAL DESIGN INC is an ISO 13485 registered company, Certificate No.: M 8072.

FlowArt is a registered trademark of Asset Medikal.

Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions.

©2021 Asset Medikal. All rights reserved.

Asset Medikal Tasarım Sanayi ve Ticaret A.Ş.

İkitelli OSB Mahallesi
17. Cadde, No: 17
Basaksehir, Istanbul - Turkey

Phone : +90 212 494 27 27

Fax : +90 212 494 34 37

Taxation Authority / No : İKİTELLİ V.D. / 0910501228

Trade Register No : 695458

Central Registration System No : 0091050122800019

info@assetmedikal.com | www.assetmedikal.com

FlowArt®

asset®
Design the future