

DESCRIPTION

FlowArt® BM01 is a connector with an integrated silicone seal valve.

FEATURES

- ☑ Microbial barrier swabable silicon stem with a split septum, flat seal to allow enhanced disinfection and drying time.
- ☑ Clear transparent housing to allow clear view into the fluid pathway for enhanced flushing.
- ☑ Compatible with TPN solutions, blood and oncology drugs.
- ☑ Suitable for high pressure infusions (up to 325 psi).
- ☑ Allows multiple accesses.
- ☑ ISO 80369-7 compatible.
- ☑ Photoprotective (-PH) versions are in amber color.
- ☑ Latex-free, PVC-free, DEHP-free

Part	Material
Valve Body	Copolyester
Swabable septum	Silicone

DIAGRAM**PRODUCT DATA**

Characteristic	Specification
Product Code	BM01, BM01-PH, BM01-M*
Product Name	Needle-free Valve Port
Priming Volume (mL)	0.085
Displacement (mL)	-0.02 (negative)
Flow Rate (mL/min @ 1m water head pressure)	312
Flow Rate (mL/min @ 0.5 bar pressure)	645
Flushing Volume (mL)	2
Length (cm)	2.5

*blue lower part

LABELLING AND PACKAGING

Units are double bagged with polyethylene bags, placed in corrugated shipper, labelled on its bag and packaging container.

Each container is labelled with the following traceability information:

- Quantity,
- Product description,
- Product date,
- Lot number

TDS.107_Rev03_26.04.2021

ASSET® MEDICAL DESIGN INC is an ISO 13485 registered company, Certificate No.: M 8072.

FlowArt is a registered trademark of Asset Medikal.

©2021 Asset Medikal. All rights reserved.

Asset Medikal Tasarım Sanayi ve Ticaret A.Ş.

İkitelli OSB Mahallesi
17. Cadde, No: 17
Basaksehir, İstanbul - Turkey

Phone : +90 212 494 27 27

Fax : +90 212 494 34 37

Taxation Authority / No : İKİTELLİ V.D. / 0910501228

Trade Register No : 695458

Central Registration System No : 0091050122800019

info@assetmedikal.com | www.assetmedikal.com

FlowArt®

asset®
Design the future