






DESCRIPTION

Usage:
on/off control of fluid or air

FEATURES

- Open Jaw
- For tubings with OD 2-3 mm
- Adjustment controls flow from full open to complete closure
- Will not distort, rust
- Materials: Polyoxymethylene (POM)

DIAGRAM**PRODUCT DATA**

Characteristic	Specification			
Product Code	K009-K	K009-M	K009-B	K009-Y
Drawing				
Color	Red	Blue	White	Green
Width (mm)	3.5			
Height (mm)	6			
Length (mm)	21			

LABELLING AND PACKAGING

Each unit is labelled on its packaging container.

Each container is labelled with the following traceability information:

- Quantity,
- Product description,
- Production date,
- Lot number

TDS.119_Rev02_26.04.2021

ASSET® MEDICAL DESIGN INC is an ISO 13485 registered company, Certificate No.: M 8072.
©2021 Asset Medikal. All rights reserved.

Asset Medikal Tasarım Sanayi ve Ticaret A.Ş.

İkitelli OSB Mahallesi
17. Cadde, No: 17
Basakşehir, İstanbul - Turkey

Phone : +90 212 494 27 27

Fax : +90 212 494 34 37

Taxation Authority / No : İKİTELLİ V.D. / 0910501228

Trade Register No : 695458

Central Registration System No : 0091050122800019

info@assetmedikal.com | www.assetmedikal.com